

# Emperor

# EmFast

Tubular Bandage

Safe and easy to use, these non-sterile, viscose elastic tubular bandages offer 2-way stretch to ensure conformability and mobility for the patients.

EmFast bandages can be used for dressing fixation, under wet or dry patch wrapping, as a skin covering or to hold ointments in place.



Dressing Fixation



Ointment Retention



Undercast Stockinette



Suitable for Wet & Dry Wrapping



Two-Way Stretch



Multiple Range of Sizes



Reusable after Washing



Can Be Cut to Required Length

For **FREE** product samples and more information, please get in touch

**Call:** 0800 008 6646

**E-mail:** [woundcare@rociallehealthcare.com](mailto:woundcare@rociallehealthcare.com)

**Visit:** [www.emperorwound.com](http://www.emperorwound.com)

 Rociale  
Healthcare

Emperor is a product brand of Rociale Healthcare



## Delivering continuity of care... ...providing value in practice!

10.75cm x 5m Yellow Line

<b>EmFast</b>	<b>£5.91</b>
Comfast	£6.04
Acti-Fast 2-Way stretch	£6.05
Tubifast 2-Way stretch	£6.09
<b>CliniFast</b>	<b>£6.23</b>

Based on DT Prices August 2024



\*Based on DT cost comparison of EmFast vs Clinifast 10.75cm x 1m Yellow Line, August 2024

### EmFast Range

Colour Descriptor	Size	Limb Circumference	Product Code	NHSSC Code	PIP Code	Drug Tariff Price
Red Line (Small)	3.5cm x 1m	9-18cm	538018ZD	EGP85014	425-1377	£0.51
	3.5 x 10m	9-18cm	538023ZD	EGP85021	-	-
Green Line (Small/Medium)	5cm x 1m	14-24cm	538020ZD	EGP85010	425-1633	£0.52
	5cm x 5m	14-24cm	538040ZD	EGP85007	425-1666	£2.73
	5cm x 10m	14-24cm	538024ZD	EGP85006	-	-
Blue Line (Large)	7.5cm x 1m	24-40cm	538021ZD	EGP85008	425-1658	£0.70
	7.5cm x 3m	24-40cm	538036ZD	EGP85016	425-1419	£2.02
	7.5cm x 5m	24-40cm	538034ZD	EGP85019	425-1641	£3.54
	7.5cm x 10m	24-40cm	538030ZD	EGP85018	-	-
Yellow Line (Extra Large)	10.75cm x 1m	35-64cm	538013ZD	EGP85008	425-1427	£1.09
	10.75cm x 3m	35-64cm	538038ZD	EGP85016	425-1385	£3.34
	10.75cm x 5m	35-64cm	538035ZD	EGP85019	425-1625	£5.91
	10.75cm x 10m	35-64cm	538031ZD	EGP85018	-	-
Beige Line (Adult)	17.5cm x 1m	64-100cm	538015ZD	EGP85012	425-1617	£1.66
	17.5cm x 10m	64-100cm	538032ZD	EGP235	-	-
Purple Line (Large Adult Trunk)	20cm x 1m	100-130cm	538017ZD	EGP85011	425-1674	£3.16
	20cm x 5m	100-130cm	538043ZD	EGP85017	425-1708	£16.09
	20cm x 10m	100-130cm	538033ZD	EGP211	-	-

**Call:** 0800 008 6646

**E-mail:** [woundcare@rociallehealthcare.com](mailto:woundcare@rociallehealthcare.com)

**Visit:** [www.emperorwound.com](http://www.emperorwound.com)



Emperor is a product brand of Rociale Healthcare