



# Emperor

## EmFilm

Film Dressing

Made from latex-free, semi-permeable polyurethane film membrane and coated with a hypoallergenic acrylic adhesive, EmFilm Dressing are designed to be breathable with good oxygen and moisture vapour exchange.

Waterproof and impermeable to liquids, bacteria and viruses, the risk of infection is reduced while also enabling effective wound monitoring through a transparent film.



Protection of newly forming tissue



Fixation of primary dressings



Provision of a 'moist wound healing environment'



The patient can shower with the dressing in place



Prevention of contamination by dirt or bacteria, thus a reduced risk of wound infection



Rounded edges and consistent hypoallergenic adhesive help prevent dressing damage and early removal, reducing the number of dressing changes needed.

**Call:** 0800 008 6646

**E-mail:** [woundcare@rociallehealthcare.com](mailto:woundcare@rociallehealthcare.com)

**Visit:** [www.emperorwound.com](http://www.emperorwound.com)




**Rociale**  
Healthcare

Emperor is a product brand of Rociale Healthcare

# Save up to 80% vs other brands\*

- EmFilm 10x15cm
- Hydrofilm 10x15cm
- ClearFilm 12x12cm
- Tegaderm 12x12cm



\*Based on DT prices August 2024

Size	Product Code	NHSSC Code	PIP Code	Drug Tariff Price
4cm x 4.4cm	537695ZD	ELW85022	425-0007	-
5cm x 7cm	537687ZD	-	425-0015	£0.10
8cm x 10cm	537696ZD	ELW85030	424-9983	£0.16
10cm x 15cm	537688ZD	-	424-9991	£0.23
10cm x 25cm	537689ZD	-	426-0311	-
20cm x 30cm	537697ZD	ELW85032	425-0056	£0.87
5cm x 10m (Non-Sterile roll)	537699ZD	ELW85024	-	-
10cm x 10m (Non-Sterile roll)	537700ZD	ELW1065	-	-

To view our full portfolio of  
Wound Care products, visit

[www.emperorwound.com](http://www.emperorwound.com)



For FREE product samples and more information, please get in touch

Call: 0800 008 6646

E-mail: [woundcare@rociallehealthcare.com](mailto:woundcare@rociallehealthcare.com)

Visit: [www.emperorwound.com](http://www.emperorwound.com)



Emperor is a product brand of Rocialle Healthcare