



# Emperor

## EmPore

Adhesive Island Dressing

Made from latex-free, porous non-woven fabric and coated with a hypoallergenic acrylic adhesive, EmPore Dressings are designed to be breathable with consistent oxygen and moisture vapour exchange.

EmPore dressing provides a protective barrier against further trauma and contamination, while its rounded corners prevent lifting decreasing the need for frequent dressing changes.



Water repellent backing, allowing the patient to shower with the dressing in place



Prevention of contamination by dirt or bacteria, thus a reduced risk of wound infection



Absorbent wound pad, designed to absorb blood and exudate



Conformable to all body parts, allowing freedom of movement



The non-woven fabric backing enables air and water vapour permeability due to its porous structure, while also conforming to the body and allowing freedom of movement to ensure patient comfort.

**Call:** 0800 008 6646

**E-mail:** [woundcare@rociallehealthcare.com](mailto:woundcare@rociallehealthcare.com)


**Visit:** [www.emperorwound.com](http://www.emperorwound.com)



**Rociale**  
Healthcare

Emperor is a product brand of Rociale Healthcare

# Supporting positive changes... ...With simple switches



Based on DT cost comparison vs comparative sizes October 2024\*

**For FREE product samples and more information, please get in touch**

Size	Product Code	NHSSC Code	PIP Code	Drug Tariff Price
5cm x 7cm	537833ZD	EIJ85012	424-3457	£0.03
8cm x 10cm	537834ZD	EIJ85009	424-3481	£0.05
10cm x 15cm	537835ZD	EIJ85008	424-3473	£0.09
10cm x 20cm	537836ZD	EIJ85010	424-3465	£0.12
10cm x 25cm	537837ZD	EIJ85007	424-3499	£0.12
10cm x 30cm	537838ZD	EIJ85011	424-3549	£0.13
10cm x 35cm	537839ZD	EIJ85006	424-2400	£0.30

**Call:** 0800 008 6646

**E-mail:** [woundcare@rocialehealthcare.com](mailto:woundcare@rocialehealthcare.com)

**Visit:** [www.emperorwound.com](http://www.emperorwound.com)



Emperor is a product brand of Rociale Healthcare